

Hematology Survey Calibration Verification /Linearity Kit




PS by Panudda Sothanapaisan

CL Rihs

20 Nov 2024

Materials For LN9



MATERIALS FOR THIS MAILING

LN9 LN9-08R — LN9-13R (RBC/HGB Set)
LN9-08W — LN9-14W (WBC Set)
LN9-08P — LN9-14P (PLT Set)



Storage and Stability

Unopened Storage

Store at 2-8°C, do not freeze.

Open Storage

Use immediately, then discard.

Appearance

Normal appearance indicates product integrity.

Visible hemolysis of the supernatant may indicate product deterioration. **Do not use product if deterioration is suspected.**



Appearance



Red blood cell/hemoglobin (RBC/HGB)

Fresh whole blood , A light pink-tinted supernatant



White blood cell (WBC)
clear to pale yellow supernatant



Platelet (PLT)
Clear to brownish supernatant

LN9-B 2024

Hematology Survey

Calibration Verification/Linearity
Kit Instructions ©CAP 2024

Hematology Calibration Verification/Linearity	
Analyte	Program Code
	LN9
Hemoglobin	■
Platelet count	■
RBC count	■
WBC count	■

Program Information

- Twenty 3.0-mL liquid specimens
- Conventional and International System of Units (SI) reporting offered

Shipping Schedule

- Shipment A: April 4
- Shipment B: September 26



Detailed Testing Instructions

RBC/HGB Set

RBC/HGB Set

- a. Allow vials to warm at room temperature for at least 15 minutes before mixing.
- b. To mix, hold vial horizontally between the palms of the hands.
- c. Roll the vial back and forth **rapidly for 30 seconds**. Gently invert the vial 10 times.
- d. Continue to mix the vial in this manner until the cells are completely suspended and there are no visible aggregates. **Do not shake the vial or mix on a mechanical (vortex) blood mixer.**

WBC/PLT Sets

WBC/PLT Sets

- a. Allow vials to warm at room temperature for at least 15 minutes before mixing.
- b. Mix by vigorously agitating on a vortex mixer for 2 minutes to ensure complete dispersion of microaggregates.
- c. Let vials sit undisturbed for 10 minutes to allow micro bubbles to dissipate before sampling.
- d. **For WBC and PLT testing on instruments with a blood detector:** Turn off the blood detector between each sample analysis. These instruments may not be able to test the specimens that do not contain any red blood cells.

1

2

3

4



Gently invert the vial 10 times immediately before each sampling.

Perform **2** assays from each vial within the same run. Two data points **must** be received for each solution.

At least **4** consecutive data sets **must** be submitted to receive an evaluation for both calibration verification and linearity.

How to Mix the Samples

RBC/HGB Set



- Allow vials to warm at room temperature for at least 15 minutes before mixing.
- Roll the vial back and forth **rapidly for 30 seconds**. Gently invert the vial 10 times.
- Continue to mix the vial in this manner until the cells are completely suspended and there are no visible aggregates.

*****Do not shake the vial or mix on a mechanical (vortex) blood mixer*****

WBC/PLT Sets

How to Mix the Samples





- Allow vials to warm at room temperature for at least 15 minutes before mixing
- Mix by vigorously agitating on a vortex mixer for 2 minutes to ensure complete dispersion of microaggregates. Let vials sit undisturbed for 10 minutes to allow micro bubbles to dissipate before sampling.

*****For WBC and PLT testing on instruments with a blood detector: Turn off the blood detector between each sample analysis.** These instruments may not be able to test the specimens that do not contain any red blood cells.

How to perform the Samples

- Gently invert the vial 10 times immediately before each sampling.
- Perform 2 assays from each vial within the same run. Two data points must be received for each solution




****For WBC and PLT testing on instruments with a blood detector:
Turn off the blood detector between each sample analysis.****

Open Cap


Press start
(Blue button)

How to Submit the Results


www.cap.org



Website



Login



Select

A screenshot of the CAP Hematology Survey form. The form is titled 'Hematology Survey (Calibration Verification/Linearity) Result Form'. It includes a 'Page 1' header with the CAP logo and 'LN9-B 2016'. The form contains a 'Report' section with a date of 'October 15, 2024'. Below this is a 'Results - Red Blood Cells - 10E12/L or 10E6/μL' section. The form is divided into two main columns: 'Left Axis' and 'Right Axis'. Each column contains a table with rows for different hematology parameters (HGB, HCT, HbA, HbF, HbS, HbT, HbE, HbF, HbS, HbT, HbE) and columns for 'Lab Value' and 'Reference Range'. The form also includes a 'Notes' section at the bottom.

Data Entry

Instrument-Specific Information

INSTRUMENT-SPECIFIC INFORMATION

Instrument/method-specific notes are provided by individual manufacturer(s). If your manufacturer did not provide the CAP with any special instructions, process specimens as you would a patient specimen or contact your manufacturer directly with any questions.

1. **Abbott Cell-Dyn 4000 and Sapphire users:** Instruments may report WBCs as nucleated red blood cells (nRBC). Users should report the sum of WBC and nRBC on the result form.
2. **Abbott Cell-Dyn 3500/3700 users:** The WBC set and the platelet set do not contain RBCs. The RBCs are used to trigger sampling volume in the instrument. The operator must allow the full instrument sample time by not removing the sample until the instrument sampling needle wash block moves down.
3. **Siemens ADVIA and Sysmex XE-series users:** Analyzers have the capability of enumerating nRBC and correcting the WBC count, report the **uncorrected** WBC results.
4. **Sysmex XN-series users:** Report the Total Nucleated Count (TNC-N) as the WBC count. To obtain the TNC-N:
 - a. Select the result from the Explorer browser screen, select the "Service Tab."
 - b. Select "WNR" from the menu to the left. Find the "TNC-N" count.
 - c. Report the "TNC-N" count for the WBC.

Report the TNC-N count for the WBC

The screenshot displays the Sysmex XN-series instrument's service menu. The top navigation bar includes icons for Menu, QC File, Work List, Rule, Explorer, and Browser. The 'Data Browser' section shows 'Modify' and 'Validate' buttons. The main display area is divided into several sections:

- Positive** (Diff, Morph, Count) and **Validated** status.
- Rule Result** (None).
- Service Data** section with a **Sampling Data** table and **Scattergram Sensitivity** table.
- Reference Data** section with **WBC-N** and **TNC-N** values.
- Cell 1** and **Laser Current** (LD driver) values.

The **Sampling Data** table shows the following values:

WNR	
917	932
927	0
904	0
937	0
954	0
943	0
966	0
930	0
8413 (*6)	

The **Scattergram Sensitivity** table shows the following values:

Parameter	Value	Unit	Parameter	Value
WNR-X	232.4	ch	WNR-WX	116
WNR-Y	216.9	ch	WNR-WY	129
NRBC-X	228.3	ch		
NRBC-Y	204.7	ch		

The **Reference Data** section shows the following values:

Parameter	Value	Unit
WBC-N	21.994	10 ³ /uL
TNC-N	53.797	10 ³ /uL
Cell 1	50482	
Laser Current	50.80	mA

The **TNC-N** value is highlighted in green, indicating it is the value to be reported for the WBC count.

Evaluation Report

College of American Pathologists 325 Waukegan Road, Northfield, Illinois 60093-2750 800-323-4040 - cap.org Laboratory Quality Solutions		CAP Number: 7080301-01 Kit #: 01 Institution: Rsch Inst for Hlth Sci Chiang Mai U Attention: Warunee Jit-aree MSc City/State: Chiang Mai, 50202		Page 1 Kit ID: 38478251 Kit Mailed: 09/23/2024 Original Evaluation: 10/18/2024 Next Mailing Date: 03/31/2025		
EVALUATION ORIGINAL		LN9-B 2024 Hematology Calibration Verification/Linearity				
Executive Summary						
Analyte		Calibration Verification		Linearity Evaluation		Page #
Red Blood Cells 10E12/L (10E6/ μ L)		Different		Linear from 0.305 to 7.500		2 - 3
Hemoglobin g/dL		Verified from 0.90 to 23.40		Linear from 0.90 to 23.40		4 - 5
White Blood Cells 10E9/L (10E3/ μ L)		Verified from 0.60 to 375.95		Linear from 0.60 to 375.95		6 - 7
Platelets 10E9/L (10E3/ μ L)		Verified from 8.0 to 3823.0		Linear from 8.0 to 3823.0		8 - 9
Troubleshooting resources: - Calibration Verification Troubleshooting Guide and Investigation Checklist - available in the Participant Summary and ELSS						

Evaluation Report

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 2

EVALUATION
ORIGINAL


LN9-B 2024 Hematology Calibration Verification/Linearity
Red Blood Cells 10E12/L (10E6/ μ L) Calibration Verification Evaluation

Evaluation Result: Different
Peer Instrument: SYSMEX XN-SERIES

Allowable Error: 3% or 0.2 10E12/L (10E6/ μ L),
whichever is greater

Specimen	Assay 1	Assay 2	Your Mean	Peer Mean	Peer N	Difference	Allowable Error
LN9-08R	0.31	0.30	0.305	0.318	571	-0.013x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-09R	1.82	1.80	1.810	1.916	574	-0.106x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-10R	3.16	3.15	3.155	3.450	574	-0.295x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-11R	4.56	4.55	4.555	5.073	574	-0.518x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-12R	5.72	5.78	5.750	6.566	574	-0.816x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-13R	7.38	7.62	7.500	8.138	574	-7.8%	$\pm 3.0\%$

Note: 10E12/L is equivalent to 10E6/ μ L



RBC

Peer Results Summary Table:
Evaluation of Instrument Performance

Peer Group Size: 574

Range	Calibration Verification		Linearity Evaluation		
	% Verified	% Different	% Linear	% Nonlinear	% Imprecise
LN9-08R - 13R	95.1	1.7	97.2	0.2	1.2
LN9-08R - 12R	1.0	0.0	0.3	0.0	0.0
LN9-09R - 13R	0.7	0.2	0.3	0.0	0.2
LN9-08R - 11R	1.2	0.0	0.5	0.0	0.0

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 4


EVALUATION
ORIGINAL

LN9-B 2024 Hematology Calibration Verification/Linearity
Hemoglobin g/dL Calibration Verification Evaluation

Evaluation Result: Verified from 0.90 to 23.40
Peer Instrument: SYSMEX XN-SERIES

Allowable Error: 3.5% or 0.67 g/dL,
whichever is greater

Specimen	Assay 1	Assay 2	Your Mean	Peer Mean	Peer N	Difference	Allowable Error
LN9-08R	0.9	0.9	0.90	0.88	567	0.02 g/dL	± 0.67 g/dL
LN9-09R	5.5	5.5	5.50	5.51	573	-0.01 g/dL	± 0.67 g/dL
LN9-10R	9.9	9.9	9.90	9.91	573	-0.01 g/dL	± 0.67 g/dL
LN9-11R	14.3	14.3	14.30	14.51	573	-0.21 g/dL	± 0.67 g/dL
LN9-12R	18.7	18.8	18.75	18.71	573	0.04 g/dL	± 0.67 g/dL
LN9-13R	23.5	23.3	23.40	23.20	572	0.9%	$\pm 3.5\%$



HGB

Peer Results Summary Table:
Evaluation of Instrument Performance

Peer Group Size: 573

Range	Calibration Verification		Linearity Evaluation		
	% Verified	% Different	% Linear	% Nonlinear	% Imprecise
LN9-08R - 13R	96.9	1.0	97.9	0.0	0.7
LN9-08R - 12R	0.9	0.0	0.3	0.0	0.0
LN9-09R - 13R	1.0	0.0	0.9	0.0	0.0
LN9-09R - 12R	0.2	0.0	0.2	0.0	0.0

Evaluation Report

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 6

EVALUATION ORIGINAL

LN9-B 2024 Hematology Calibration Verification/Linearity White Blood Cells 10E9/L (10E3/ μ L) Calibration Verification Evaluation

Evaluation Result: Verified from 0.60 to 375.95


Peer Instrument: SYSMEX XN-SERIES

Allowable Error: 7.5% or 0.12 10E9/L (10E3/ μ L),
whichever is greater

Specimen	Assay 1	Assay 2	Your Mean	Peer Mean	Peer N	Difference	Allowable Error
LN9-08W	0.6	0.6	0.60	0.60	568	0.00x10E9/L	$\pm 0.12 \times 10E9/L$
LN9-09W	10.9	10.9	10.90	10.58	569	3.0%	$\pm 7.5\%$
LN9-10W	53.8	52.9	53.35	52.10	569	2.4%	$\pm 7.5\%$
LN9-11W	100.0	100.7	100.35	97.70	570	2.7%	$\pm 7.5\%$
LN9-12W	198.8	198.1	198.45	196.95	569	0.8%	$\pm 7.5\%$
LN9-13W	305.3	303.1	304.20	298.65	569	1.9%	$\pm 7.5\%$
LN9-14W	374.8	377.1	375.95	368.36	569	2.1%	$\pm 7.5\%$

Note: 10E9/L is equivalent to 10E3/ μ L

Calibration Verification Plot: Percent Differences with Allowable Error Limits



WBC

Peer Results Summary Table:
Evaluation of Instrument Performance

Peer Group Size: 570

Range	Calibration Verification		Linearity Evaluation		
	% Verified	% Different	% Linear	% Nonlinear	% Imprecise
LN9-08W - 14W	80.9	14.7	87.9	0.4	2.1
LN9-08W - 13W	0.5	0.0	1.1	0.0	0.0
LN9-09W - 14W	0.9	0.0	0.2	0.0	0.0
LN9-08W - 12W	0.4	0.0	2.3	0.0	0.0
LN9-08W - 11W	2.3	0.2	6.0	0.0	0.0
LN9-11W - 14W	0.2	0.0	0.2	0.0	0.0

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 8

EVALUATION ORIGINAL

LN9-B 2024 Hematology Calibration Verification/Linearity Platelets 10E9/L (10E3/ μ L) Calibration Verification Evaluation

Evaluation Result: Verified from 8.0 to 3823.0


Peer Instrument: SYSMEX XN-SERIES

Allowable Error: 12.5% or 4.5 10E9/L (10E3/ μ L),
whichever is greater

Specimen	Assay 1	Assay 2	Your Mean	Peer Mean	Peer N	Difference	Allowable Error
LN9-08P	8	8	8.0	7.9	572	0.1x10E9/L	$\pm 4.5 \times 10E9/L$
LN9-09P	46	45	45.5	51.1	573	-11.0%	$\pm 12.5\%$
LN9-10P	371	366	368.5	391.7	573	-5.9%	$\pm 12.5\%$
LN9-11P	741	729	735.0	786.7	573	-6.8%	$\pm 12.5\%$
LN9-12P	1029	1028	1028.5	1113.4	570	-7.8%	$\pm 12.5\%$
LN9-13P	2073	2105	2089.0	2241.2	570	-6.8%	$\pm 12.5\%$
LN9-14P	3833	3813	3823.0	4116.5	569	-7.1%	$\pm 12.5\%$

Note: 10E9/L is equivalent to 10E3/ μ L

Calibration Verification Plot: Percent Differences with Allowable Error Limits



PLT

Peer Results Summary Table:
Evaluation of Instrument Performance

Peer Group Size: 573

Range	Calibration Verification		Linearity Evaluation		
	% Verified	% Different	% Linear	% Nonlinear	% Imprecise
LN9-08P - 14P	91.6	5.6	97.9	0.2	0.0
LN9-08P - 13P	1.0	0.2	1.0	0.0	0.0
LN9-09P - 14P	0.7	0.0	0.2	0.0	0.0
LN9-08P - 12P	0.2	0.0	0.2	0.0	0.0
LN9-08P - 11P	0.7	0.0	0.5	0.0	0.0

Evaluation Report

Page 3

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 • cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01 Kit ID: 38478251
Institution: Rsch Inst for Hlth Sci Chiang Mai U Kit Mailed: 09/23/2024
Attention: Warunee Jit-aree MSc Original Evaluation: 10/18/2024
City/State: Chiang Mai, 50202 Next Mailing Date: 03/31/2025

EVALUATION ORIGINAL

LN9-B 2024 Hematology Calibration Verification/Linearity Red Blood Cells 10E12/L (10E6/ μ L) Linearity Evaluation


Evaluation Result: Linear from 0.305 to 7.500

Instrument: SYSMEX XN-SERIES

Evaluation Type: Standard
Goal for Total Error (TE): 6%


Specimen	Assay 1	Assay 2	Your Mean	Best-fit Target	Relative Concentration
LN9-08R	0.31	0.30	0.305	0.295	0.000
LN9-09R	1.82	1.80	1.810	1.765	0.208
LN9-10R	3.16	3.15	3.155	3.144	0.403
LN9-11R	4.56	4.55	4.555	4.614	0.611
LN9-12R	5.72	5.78	5.750	5.894	0.792
LN9-13R	7.38	7.62	7.500	7.384	1.000

Linearity Plot 1: Reported Results with Best-fit Line



RBC

Linearity Plot 2: Differences with Limits of Acceptable Imprecision



Mean of Included Results: 3.846 10E12/L (10E6/ μ L)

◆ Included in best-fit line
○ Excluded from best-fit line

Your plot has one or more points within your linear range that fall outside of the shaded area. Since your evaluation is Linear, no remedial action is necessary.

Points can fall outside of the shaded area for two reasons:
1) an average is used to estimate imprecision, so many small differences can offset a few large differences, and
2) clinically insignificant nonlinearity (curved fit) can contribute to differences between your results and the best-fit straight line. Larger differences may be an early warning sign of nonlinearity, poor repeatability, or poor fit.

Page 5

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 • cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01 Kit ID: 38478251
Institution: Rsch Inst for Hlth Sci Chiang Mai U Kit Mailed: 09/23/2024
Attention: Warunee Jit-aree MSc Original Evaluation: 10/18/2024
City/State: Chiang Mai, 50202 Next Mailing Date: 03/31/2025

EVALUATION ORIGINAL

LN9-B 2024 Hematology Calibration Verification/Linearity Hemoglobin g/dL Linearity Evaluation


Evaluation Result: Linear from 0.90 to 23.40

Instrument: SYSMEX XN-SERIES

Evaluation Type: Standard
Goal for Total Error (TE): 7%


Specimen	Assay 1	Assay 2	Your Mean	Best-fit Target	Relative Concentration
LN9-08R	0.9	0.9	0.90	0.82	0.000
LN9-09R	5.5	5.5	5.50	5.50	0.208
LN9-10R	9.9	9.9	9.90	9.89	0.403
LN9-11R	14.3	14.3	14.30	14.57	0.611
LN9-12R	18.7	18.8	18.75	18.64	0.792
LN9-13R	23.5	23.3	23.40	23.32	1.000

Linearity Plot 1: Reported Results with Best-fit Line



HGB

Linearity Plot 2: Differences with Limits of Acceptable Imprecision



Mean of Included Results: 12.12 g/dL

◆ Included in best-fit line
○ Excluded from best-fit line

Evaluation Report

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 7

EVALUATION
ORIGINAL


LN9-B 2024 Hematology Calibration Verification/Linearity
White Blood Cells 10E9/L (10E3/ μ L) Linearity Evaluation

Evaluation Result: Linear from 0.60 to 375.95
Instrument: SYSMEX XN-SERIES


Evaluation Type: Standard
Goal for Total Error (TE): 15%

Specimen	Assay 1	Assay 2	Your Mean	Best-fit Target	Relative Concentration
LN9-08W	0.6	0.6	0.60	4.41	0.000
LN9-09W	10.9	10.9	10.90	13.39	0.024
LN9-10W	53.8	52.9	53.35	51.93	0.127
LN9-11W	100.0	100.7	100.35	94.96	0.242
LN9-12W	198.8	198.1	198.45	196.74	0.514
LN9-13W	305.3	303.1	304.20	303.76	0.800
LN9-14W	374.8	377.1	375.95	378.60	1.000

Linearity Plot 1: Reported Results with Best-fit Line



Linearity Plot 2: Differences with Limits of Acceptable Imprecision



Mean of Included Results: 149.11 10E9/L (10E3/ μ L)

WBC

◆ Included in best-fit line
○ Excluded from best-fit line

College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 - cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01
Institution: Rsch Inst for Hlth Sci Chiang Mai U
Attention: Warunee Jit-aree MSc
City/State: Chiang Mai, 50202

Kit ID: 38478251
Kit Mailed: 09/23/2024
Original Evaluation: 10/18/2024
Next Mailing Date: 03/31/2025

Page 9

EVALUATION
ORIGINAL


LN9-B 2024 Hematology Calibration Verification/Linearity
Platelets 10E9/L (10E3/ μ L) Linearity Evaluation

Evaluation Result: Linear from 8.0 to 3823.0
Instrument: SYSMEX XN-SERIES


Evaluation Type: Standard
Goal for Total Error (TE): 25%

Specimen	Assay 1	Assay 2	Your Mean	Best-fit Target	Relative Concentration
LN9-08P	8	8	8.0	< 0	0.000
LN9-09P	46	45	45.5	32.8	0.011
LN9-10P	371	366	368.5	364.6	0.098
LN9-11P	741	729	735.0	746.0	0.198
LN9-12P	1029	1028	1028.5	1051.1	0.278
LN9-13P	2073	2105	2089.0	2107.6	0.555
LN9-14P	3833	3813	3823.0	3804.8	1.000

Linearity Plot 1: Reported Results with Best-fit Line



Linearity Plot 2: Differences with Limits of Acceptable Imprecision



Mean of Included Results: 1156.8 10E9/L (10E3/ μ L)

PLT

◆ Included in best-fit line
○ Excluded from best-fit line

Troubleshooting

Page 2
College of American Pathologists
325 Waukegan Road, Northfield, Illinois 60093-2750
800-323-4040 • cap.org
Laboratory Quality Solutions

CAP Number: 7080301-01 Kit #: 01 Kit ID: 38478251
Institution: Rsch Inst for Hlth Sci Chiang Mai U Kit Mailed: 09/23/2024
Attention: Warunee Jit-aree MSc Original Evaluation: 10/18/2024
City/State: Chiang Mai, 50202 Next Mailing Date: 03/31/2025

EVALUATION
ORIGINAL


LN9-B 2024 Hematology Calibration Verification/Linearity
Red Blood Cells 10E12/L (10E6/ μ L) Calibration Verification Evaluation

Evaluation Result: Different
Peer Instrument: SYSMEX XN-SERIES

Allowable Error: 3% or 0.2 10E12/L (10E6/ μ L),
whichever is greater

Specimen	Assay 1	Assay 2	Your Mean	Peer Mean	Peer N	Difference	Allowable Error
LN9-08R	0.31	0.30	0.305	0.318	571	-0.013x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-09R	1.82	1.80	1.810	1.916	574	-0.106x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-10R	3.16	3.15	3.155	3.450	574	-0.295x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-11R	4.56	4.55	4.555	5.073	574	-0.518x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-12R	5.72	5.78	5.750	6.566	574	-0.816x10E12/L	$\pm 0.200 \times 10E12/L$
LN9-13R	7.38	7.62	7.500	8.138	574	-7.8%	$\pm 3.0\%$

Note: 10E12/L is equivalent to 10E6/ μ L.



Peer Results Summary Table:
Evaluation of Instrument Performance


Peer Group Size: 574

Range	Calibration Verification		Linearity Evaluation		
	% Verified	% Different	% Linear	% Nonlinear	% Imprecise
LN9-08R - 13R	95.1	1.7	97.2	0.2	1.2
LN9-08R - 12R	1.0	0.0	0.3	0.0	0.0
LN9-09R - 13R	0.7	0.2	0.3	0.0	0.2
LN9-08R - 11R	1.2	0.0	0.5	0.0	0.0

Calibration Verification Troubleshooting Guide

This troubleshooting guide provides suggested actions if you receive a calibration verification evaluation result of Different, or if your evaluation result is Verified over a range that does not include all of your reported results. To use this guide, determine which of the following examples is most similar to your calibration verification plot. Refer to the corresponding suggested actions, in conjunction with the CVL Investigation Checklist for Problematic Results, to investigate possible causes and corrective actions.

Constant Bias



Review the **ANALYTICAL** section of the investigation checklist. Analytical problems that produce a constant bias may be due to a calibration error. Recalibration may be needed.

Review the **CLERICAL** section of the investigation checklist. Clerical errors that result in a constant bias are likely due to units of measure or decimal place errors, or incorrect peer group assignment. Clerical errors may indicate a need for additional staff training.

CVL Survey Investigation Checklist for Problematic Results

Analytical

YES NO NA

Was the written procedure followed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was instrument maintenance performed on schedule?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Were quality control results acceptable?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the most recent calibration acceptable and within established stability limits at the time testing was performed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Does a review of recent proficiency testing results or past CVL results indicate evenly distributed data without bias?	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Were the reagents prepared according to procedure?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Were the reagents within their open stability acceptable range?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the intended result within the measuring range for the instrument?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the dilution protocol followed when diluting samples that are out of range?	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does a review of records indicate that there were no related instrument/test problems noted prior to or after the testing was performed?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

A response of "No" to any of these questions may indicate an analytical error. These types of errors could indicate a failure to follow recommended instrument maintenance and calibration. You may need to review the instructions provided with the testing material and/or laboratory procedures. If recalibration has not already occurred, recalibrate the instrument.

CVL Survey Investigation Checklist for Problematic Results

Clerical

YES NO NA

Were the results correctly transcribed from the instrument read-out or report?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Was the correct instrument/method/reagent code reported on the result form?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Do the units of measure match between the result form and the instrument results?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Is the decimal place correct?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Does the submitted result match the result found on the calibration verification evaluation report?	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
If the result was out of range and a dilution was performed, was the correct dilution factor used in the calculation of the final result?	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

A response of "No" to any of these questions may indicate a clerical error. Although reporting of testing results is unlike those for patient results, clerical errors may indicate a need for additional staff training, review of kit instructions, or investigation of the reporting format provided by the testing device. If results reported on the result form do not match the results found on the evaluation report, please contact the CAP Customer Contact Center at 800-323-4040.

EQA Evaluation for Quantitative Testing

Date:	28 Aug 2024	EQA Provider:	OWA: TH1021AI
Laboratory:	HID0131_RIHES Research Institute for Health Sciences at Chiang Mai University, Chiang Mai, Thailand	Panel:	Hematology HEFG435-2 2024

Summary: This EQA event was successful for all protocol analytes.

Investigation Reports are required for all **protocol analytes** scoring less than 100%. The analytes requiring investigation are:
None

Internal investigations are recommended for non-protocol analytes scoring less than 100% , and for any bias, shifts and/or trends identified below:

Bias noted:

Positive: Platelet Count
Negative: None

Shifts noted:

Positive: None
Negative: None

Trends noted:

Positive: None
Negative: None

Non-Protocol Analytes scoring less than 100%:
None

Comments:


Congratulations on your successful EQA results.

Reviewer:

Amy Rada, BS MLS(ASCP)CM
Senior International QA/QC Coordinator


W. J. Rada
Digitally signed by WJ Rada
DN: cn=WJ Rada, o=Chiang Mai
University, cn=WJ Rada, c=Thailand
Reason: I have reviewed all pages
of this document.
Location: Clinical Lab, RHES CMU
Date: 2024.09.02 10:52:11 +0700
Foxit PDF Editor Version: 2024.3.0

Performance Report


Participant	TH1021AI	Research Institute for Health Sciences
Subscription ID	275066	HEFG435 Hematology 5-Part Differential
Results Deadline	2024/Jul/17	
Accreditation		This program is provided by Oneworld Accuracy Inc. within the scope of its ISO/IEC 17043:2010 accreditation by A2LA under certificate 4839.01.

Instrument Model: Sysmex XN-1000

PARTICIPANT DATA					COMPARATIVE GROUP DATA						
Analyte / Sample	Result	PAD Score(%)	Grade	Statistical Count	Mean	SD	CV(%)	PG:Description	Acceptable Range	Evaluation Criteria	
Red Blood Cell Count 10^6/ L											
A	4.74	11.86	ACC	41	4.712	0.046	1.0	ID:Sysmex XN-1000	4.47 - 4.95	Peer Group Mean	5.0%
B	5.71	-8.71	ACC	41	5.735	0.055	1.0	ID:Sysmex XN-1000	5.44 - 6.03	Peer Group Mean	5.0%
C	4.69	-5.96	ACC	39	4.704	0.042	0.9	ID:Sysmex XN-1000	4.46 - 4.94	Peer Group Mean	5.0%
D	5.74	6.64	ACC	39	5.721	0.057	1.0	ID:Sysmex XN-1000	5.43 - 6.01	Peer Group Mean	5.0%
E	4.70	0.85	ACC	39	4.698	0.049	1.0	ID:Sysmex XN-1000	4.46 - 4.94	Peer Group Mean	5.0%



Monthly IQC September 2024



Monthly IQC October 2024


Clinical Laboratory, RIHES, CMU XN series admin

QC Chart

Shift All

Nickname	Material	Lot No.	Exp. Date	Date From	To
XN-1000-1-A	Control Level1	QC-42121101	20/10/2024	01/10/2024	18/10/2024

UL	SD
Target	Mean
LL	CV (%)



01/10/2024 09:10
02/10/2024 08:53
03/10/2024 09:52
04/10/2024 09:09
05/10/2024 09:14
06/10/2024 09:35
07/10/2024 08:56
08/10/2024 08:50
09/10/2024 08:59
10/10/2024 08:59
11/10/2024 08:59
12/10/2024 08:59
13/10/2024 08:59
14/10/2024 08:59
15/10/2024 08:59
16/10/2024 08:59
17/10/2024 08:59
18/10/2024 08:59

THANK
YOU FOR YOUR
ATTENTION!